



LoocidBCP™ Drill Protocols

LoocidBCP™ Drill Bits used for Nobel Biocare implants NobelParallel CC

NobelParallel Ø	Loocid Starter Drill Bit (Optional)	Soft Bone/ Ø Loocid	Soft Bone/ Ø Parallel	Medium Bone/ Ø Loocid	Medium Bone/ Ø Parallel	Very Hard Bone/ Ø Loocid	Very Hard Bone/ Ø Parallel
NP 3.75mm	2.0mm	2.8mm	2.8mm	2.8mm	2.8/3.75mm	3.2mm	3.2/3.75mm
RP 4.3mm	2.0mm	2.8mm	2.8(3.6)mm	3.2mm	3.6/4.3mm	3.2mm	3.6 /4.3mm
RP 5.0mm	2.0mm	3.2mm	3.6mm	4.2mm	4.2/5.0mm	4.2mm	4.2 (5.0)mm
WP 5.5mm	2.0mm	4.2mm	4.6mm	4.2mm	5.0/5.5mm	4.2mm	5.0/5.5mm

The LoocidBCP™ Drill Bits can be used in conjunction with the NobelParallel Conical Connection Implant system. The surgical sequence is:

- 1) LoocidBCP™ Starter Drill Bit:
2.0mm (optional use only to verify drilling direction)
- 2) LoocidBCP™ Drill Bit:
The size of the LoocidBCP™ Drill Bit used after the optional Loocid 2.0mm Drill Bit or as the first drill bit for the implant bed preparation is shown in the table above
- 3) Nobel Biocare NobelParallel Conical Connection Drill Bit:
The Nobel Biocare Drill Bit following the use of the LoocidBCP™ Drill Bit is shown in the table above. Please finish the implant bed preparation used for the correct implant diameter and bone density by following the regular NobelParallel Drill Bit protocol

